



Ultraflux®
AV 400S / AV 600S / AV 1000S
Steam-sterilised Haemofilter



Continuous Renal Replacement Therapy

Continuous Renal Replacement Treatment (CRRT) serves for the therapy of acute renal failure and is a well-tolerated and effective method. CRRT is used to correct hydrate-, electrolyte-, acid/base metabolism and for the removal of urea and creatinine ^(1, 2, 3, 4, 5, 6).

In comparison with discontinuous therapies like intermittent hemodialysis, continuous renal replacement therapies have the decisive advantage of greater cardiovascular stability of the patient through a constant and slow elimination of fluids and substances usually eliminated with the urine. Therefore CRRT is indicated especially for cardiovascular unstable patients. CRRT furthermore allows a generous fluid supply within the scope of the parenteral nutrition.

Depending on the clinical condition of the patient, various continuous treatment modes are available:

- CAVH, continuous arterio-venous haemofiltration
- CVVH /HV-CVVH, continuous veno-venous haemofiltration/high-volume haemofiltration⁽⁷⁾
- CVVHDF, continuous veno-venous haemodiafiltration
- CVVHD, continuous veno-venous haemodialysis
- SCUF, slow continuous ultrafiltration

The Ultraflux®-series


The Ultraflux®-filters contain a Fresenius Polysulfone® membrane specially developed for continuous renal replacement therapy. With this membrane, substances with a molecular weight of up to approximately 30,000 Dalton can be eliminated, but plasma proteins like albumin, bigger molecules and cellular blood constituents are retained. The filtration characteristics of this membrane comes close to the natural human kidney's glomeruli.

Due to their high mechanical stability and excellent filter lifetimes, the Ultraflux®-filters AV 400S, AV 600S and AV 1000S are tailor-made for use in all continuous renal replacement treatments.


The large inner lumen and thin wall thickness of the capillaries in the Ultraflux®-filters improve the flow geometry and give the following advantages:

- high fluid exchange possible
- unimpeded ultrafiltrate flow
- low flow resistance in the capillaries reduces friction and shear forces
- excellent filter lifetime

Mortality




CVVH (Postdilution)



Dosage at CVVH: latest studies show a reduced mortality with an increase of the exchange volume over 20 ml/h/kg body weight. For a patient who weighs e.h. 70 kg an exchange volume of at least 2l/h is currently recommended. (adapted after Ronco et al., 2000)

Recommended filter usage

Steam-sterilised



Fresenius Polysulfone® Filter for Continuous Renal Replacement Therapy

In the course of the INLINE-steam sterilisation, the blood- and filtrate/dialysate-compartments of the filter are continuously rinsed with 121 °C hot, sterile steam for 15 minutes and dried thereafter. This process ensures a sterile and pyrogen-free product with convincing advantages:

- residues of toxic sterilisation agents like e.g. ethylene oxide are not present and cannot affect the patient
- dry delivery (reduced weight)
- no time-consuming rinsing of the filter after filling
- cost reduction through reduced saline consumption

Haemo-compatibility

The Fresenius Polysulfone® Membrane

The excellent blood compatibility of the Fresenius Polysulfone® membrane has been demonstrated in numerous publications. The following advantages have been found next to others:

- no interaction with drugs (ACE-inhibitors)
- negligible activation of the complement system
- no increased cytokine synthesis
- low activation of the coagulation system

The use of anticoagulants can be minimised which is of particular importance for patients with coagulation disorders or postoperative respectively polytraumatised patients.




An important safety aspect is the 100 % leakage testing of all filters, possible only because of the special structure of the Fresenius Polysulfone® membrane and the absence of pore filling agents like glycerine.

All production steps from the manufacturing of the membrane up to the finished dialyser are adjusted to each other resulting in constant highest quality, proven by the CE-mark.

Performance data/Technical data

	AV 400S	AV 600S	AV 1000S
Effective surface area (m ²)	0.75	1.4	1.8
Wall thickness/inner diameter (µm)	35/220	35/220	35/220
Blood filling volume (ml)	52	100	130
Recommended blood flow (ml/min)	50 – 200	100 – 350	200 – 500
Recommended filtrate flow (postdilution)	maximum 20 % of effective blood flow		
Membrane material	Fresenius Polysulfone®		
Housing material	Polycarbonate		
Potting compound	Polyurethane		
Sterilisation method	INLINE-Steam		
Art.-No.	5007341	5007361	5008981



Literature

1. Thews G: Atemgastransport und Säure-Basen-Status des Blutes. in Schmidt RF, Thews G (Hrsg) Physiologie des Menschen, 27. Auflage, Springer, 1997
2. Barenbrock M, Hausberg M, Matzkies F, de la Motte S, Schaefer RM: Effects of bicarbonate- and lactate-buffered replacement fluids on cardiovascular outcome in CVVH patients. Kidney International, 58: 1751–1757, 2000
3. Kierdorf HP, Leue C, Arns S: Lactate- or bicarbonate-buffered solutions in continuous extracorporeal renal replacement therapies. Kidney International 56 [Suppl. 72]: 32–36, 1999
4. Heering P, Ivens K, Thümer O, Brause M, Grabensee B: Acid-base balance and substitution fluid during continuous hemofiltration. Kidney International 56 [Suppl. 72]: 37–40, 1999
5. Macias WL, Clark WR: Acid-base balance in continuous renal replacement therapy. Seminars in Dialysis 9: 145–151, 1996
6. Ronco C. et al.: Effects of different doses in continuous veno-venous haemofiltration on outcomes of acute renal failure: a prospective randomised trial, The Lancet, 356: 26–30, 2000
7. Cole L, Bellomo R, Journois D, Davenport P, Baldwin I, Tipping P: High volume hemofiltration in human septic shock. Int Care Med 27: 978–986, 2001



Fresenius Medical Care