

bibag On-line Dry Bicarbonate Concentrate



bibag – the on-line dry bicarbonate concentrate for haemodialysis

for Fresenius Medical Care haemodialysis machines

The best solution

Use all the advantages, that bibag provides:

Dry concentrate

- Small and lightweight
 - Easy handling
 - Small storage requirements
 - Easy transportation
- Chemically stable
 - Long shelf life

On-Line dry concentrate

- Dissolved directly at the machine immediately before use
 - Prevention of microbiological growth
 - Safe dialysis treatment

Unique additional advantages

- Flexible bag
 - Minimal waste volume
 - Comfortable to dispose of
 - Environmentally friendly
- Automatic Emptying Programme after treatment
 - No spilling or dripping during disconnection
 - No manual removal of remaining water
 - Most hygienic way to work



Specifications

| Packaging: | units/carton | cartons/pallet |
|-------------|--------------|----------------|
| bibag 650 g | 16 | 56 |
| bibag 900 g | 12 | 56 |

Material bibag: PP/PE, Overwrap: PE, Shelf life: 3 years, CE-marked

Order information

| | Weight/Bag | Art. No. |
|-------------|------------|-----------|
| bibag | 650 g | 508 992 1 |
| bibag | 900 g | 508 991 1 |
| bibag 5008* | 650 g | 506 078 1 |
| bibag 5008* | 900 g | 506 080 1 |

*bibag 5008 is the standard module for the 5008, 5008S and 4008S classic HD machines



Fresenius Medical Care